



Product Code	31000
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279310009
PARAF	903625440

## Description

### BIOLOGICAL MICROSCOPE - 40X to 1000X

- 30° compensation free binocular head
- eyepiece 10X (Field 18 mm)
- achromatic objectives 4X, 10X, 40X(S), 100X (S, oil)
- interpupillary distance adjustment: from 55 to 75 mm
- mechanical stage, double layers (size: 125x135 mm, move range: 75x35 mm)
- Abbe condenser with iris diaphragm and filter
- halogen lamp 6 V, 20 W - brightness adjustment
- power supply: 110-220 V, 50 Hz

Quadruple nosepiece and objectives

Fine focusing knob minimum division 0.002 mm

Manual: GB, FR, IT, ES, PT, DE, GR, Arabic