



Product Code	31002
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279310023

## Description

### LED BIOLOGICAL MICROSCOPE - 40X to 1600X

High level microscope with high resolution and transmission

- infinite optical system
- wide field eyepiece WF 10X/20
- infinite semi-plan achromatic objectives 4X, 10X, 40X, 100X
- Abbe condenser
- backward quadruple nosepiece
- double layers mechanical stage (size: 150x139 mm, move range 75x52 mm)
- adjustable brightness
- power supply: 110-220 V, 50 Hz

Focusing system: coaxial coarse and fine adjustment, fine division 0.002 mm, moving range 25 mm

Large operating space with backwards nosepiece.

Sharp Image with high quality objectives.

Fine focusing knob

Manual: GB, FR, IT, ES, PT, DE