



Product Code	31145
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2511989
CND	Z12040108
EAN/UPC	4053755201594

Description

HEINE DELTA 30 PRO DERMATOSCOPE - K-235.28.305

The repositioned LEDHQ elements allow light to penetrate deeper skin layers, providing a vivid, spatial impression of lesions. This 3D-like imaging performance establishes the DELTA 30 PRO as a high-end dermatoscope. Furthermore, HEINE colorSHIFT offers a 4-stage adjustable colour temperature, from warm white to bluish white, according to clinical preference.

- ergonomic handling due to angled product design

- 10X magnification

- 3 levels of light intensity

- large field of view with real Ø 30 mm (lens 32 mm)

- crystal clear image, no matter if you prefer a closer or greater working distance

- fast change from contact to noncontact examination

- no cross-contamination: autoclavable contact plate with scale

- toggle function: switch between polarised and non-polarised illumination at the touch of a button

- connection to free HEINE DERM App

Supplied with contact plate with scale, rechargeable Li-ion battery, E4-USBC power supply with cable and case.

Technical Specifications

Diopter: -4 to +4 dpt

Illuminance: >22.000 lux in 15 mm distance without contact plate

Colour temp: 5.300 K to 6.500 K to 8.300 K to 11.000 K

Colour rendering index (CRI): ≥ 85

Working distance: 15 mm (non contact mode) or contact to skin

LED life: > 50.000 h

Power supply: Li-ion built-in rechargeable battery , 6 W

Input: USB-C: 5 V DC, 1.2 A

Charging time: 150 min

Size - weight: 195x55xh 70 mm - 0.275 kg

Protection: IP20