



Product Code	31146
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2790695
CND	Z12040108
EAN/UPC	4053755197163

Description

HEINE DELTA 30 DERMATOSCOPE - K-230.28.305

Delta 30 dermatoscope stands out due to an excellent colour rendering in LEDHQ, a large field of view and an extremely sharp image with new achromatic optical system. Polarisation provides glare-free and nonreflective working conditions. It is comfortable to operate with its angled, ergonomic design and intuitive control.

- 10X magnification
- large field of view with real Ø 30mm (lens 32mm)
- crystal clear image, no matter if you prefer a closer or greater working distance
- 3 levels of light intensity
- fast change from contact to noncontact examination
- no cross-contamination: autoclavable contact plate with scale
- toggle function: switch between polarised and non-polarised illumination made at the touch of a button
- connection to free HEINE DERM App

Supplied with contact plate with scale, rechargeable Li-ion battery, USB power supply and case.

Technical Specifications

Diopter: -4 to +4 dpt
 Illuminance: >25.000 lux in 15 mm distance without contact plate
 Colour temp: 5.700 K
 Colour rendering index (CRI): = 90
 Working distance: 15 mm (non contact mode) or contact to skin
 LED life: > 50.000 h
 Power supply: Li-ion built-in rechargeable battery , 6 W
 Input: USB 2.0 Micro B, 5V DC, 1.2 A or charging station
 Charging time: 120 min
 Size: 195x55xh 70 mm
 Protection: IP20