

HEINE MINI 3000 LED DERMATOSCOPE with soft pouch,
contact plate with scale D-888.78.021



Product Code	31157
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

EAN/UPC	4053755179121
---------	---------------

Description

MINI 3000 LEDHQ - Black
HEINE code: D-888.78.021
Contact plate with scale

Modern pocket dermatoscopes with handy dimensions and metal head.
Up to 10 hours operation time and bright LED illumination.
Twice as bright as a conventional XHL Xenon Halogen dermatoscope with maintenance-free LEDHQ illumination - no longer a need to replace the light source.
10 x magnification with focusing optics.
Exclusive Battery Performance Indicator: light intensity is regulated down to the level that corresponds to the current charge status of the battery.
Fade-Out Feature: brightness reduces slowly even with low residual capacity, for even more efficient use of dry or rechargeable batteries 20,000 ON / OFF switch cycles guaranteed.
Distortion-free image with high resolution edge to edge thanks to high-quality, colour-neutral achromatic optical system.
Special small contact plate for the examination of difficult-to-access lesions.

Supplied with handle with batteries and contact plate with scale, zipper case.

Made in Germany.