

HEINE MINI 3000 LED DERMATOSCOPE with contact plate - black



Product Code	31160
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	Yes
EAN/UPC	4053755158829
PARAF	965528641

Description

MINI 3000 LEDHQ - Black
HEINE code: D-008.78.109
Contact plate with scale

Modern pocket dermatoscopes with handy dimensions and metal head.

Up to 10 hours operation time and bright LED illumination.

Twice as bright as a conventional XHL Xenon Halogen dermatoscope with maintenance-free LEDHQ illumination - no longer a need to replace the light source.

10 x magnification with focusing optics.

Exclusive Battery Performance Indicator: light intensity is regulated down to the level that corresponds to the current charge status of the battery.

Fade-Out Feature: brightness reduces slowly even with low residual capacity, for even more efficient use of dry or rechargeable batteries 20,000 ON / OFF switch cycles guaranteed.

Distortion-free image with high resolution edge to edge thanks to high-quality, colour-neutral achromatic optical system.

Special small contact plate for the examination of difficult-to-access lesions.

Supplied with handle with batteries and contact plate with scale.

Made in Germany.