

HEINE MINI 3000 DERMATOSCOPE



Product Code	31184
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	84805
CND	Z12040108
EAN/UPC	4053755124862
UMDNS	18-021
PARAF	924800511

Description

HEINE MINI 3000 XENON-HALOGEN DERMATOSCOPE 10X - D-001.78.107

Pocket dermatoscope with high-quality optics, 10X magnification with focusing optics, metal head.

Xenon Halogen bulb.
Concentrated, bright light for uniform illumination with exact colour rendering.
High quality colour-neutral optic system with acromatic lenses, provides sharp images.

Contact plate without scale.

Supplied with dermoscopic compendium.

Operates with AA batteries or rechargeable batteries (31693) + Mini 3000 NT charger (31692) not included.