



Product Code	31187
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2807141
CND	Z12040108
EAN/UPC	8023279311877
PARAF	903119410
GMDN	18021

## Description

GIMA xenon halogen DERMATOSCOPE 2.5 V - 10X

Advanced skin surface microscope offers long-lasting halogen illumination for the true tissue colour and 10X magnification for greater details, allowing accurate diagnoses of pigmented skin lesion.

Supplied in a bag, with spare autoclavable contact plate with 0-10 mm scale.  
Latex free.

Plug in and snap, easy locking  
Supplied with 2 autoclavable contact plates (one millimetrated) Ø 2.5 cm  
Satin finish  
Multilingual manual and box (GB, FR, IT, ES, PT, SE, HU, RO)

## Technical Specifications

Battery: 2x"C" or rechargeable battery (34467) + charger (31542)  
Reticle contact plate  
Easy to focus, with large focus ring