



| | |
|---------------|----------------|
| Product Code | 31196 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2254706 |
| CND | Z129004 |
| EAN/UPC | 8436562860035 |

Description

VISTA WOOD LAMP 9W - handle

Lamps to diagnose the condition of skin that appears in different colours under ultraviolet rays. Polished glass magnifying lens 10x6 cm - 3 diopres.

Fluorescent lamps have 9W special compact fluorescent tube at approximately 360 nm wavelength.

Voltage: 220-240 V - 50/60 Hz
Manual in GB, FR, IT, ES, PT.

Made in Spain.