

VISTA WOOD LAMP - trolley



Product Code	31197
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	1887267
CND	Z120614
EAN/UPC	8436562860028

Description

VISTA WOOD LAMP 9W - trolley

Lamps to diagnose the condition of skin that appears in different colours under ultraviolet rays. Polished glass magnifying lens 10x6 cm - 3 dioptries.

Fluorescent lamps have 9W special compact fluorescent tube at approximately 360 nm wavelength.

Voltage: 220-240 V - 50/60 Hz
Manual in GB, FR, IT, ES, PT.

Arm length 45+45 cm
Trolley height 95 cm

Made in Spain.