



Product Code	31220
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

CND	Z12120122
EAN/UPC	8023279312201
PARAF	903129474

Description

SCHIÖTZ TONOMETER - straight scale

A masterpiece of precision mechanism, measures intraocular pressure according to Prof. Schiötz standards.

Vital parts are made of stainless steel.

High quality agate bearing for extremely long service life.

Set 3 weights included (5.5-7.5-10 g).

Scale specification 5.

In black deluxe case with velvet-look inserts.

Perfect reading of the scale with red pointer.

Made in Germany.