


SCHIÖTZ TONOMETER - inclined



Product Code	31222
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

CND	Z12120122
EAN/UPC	8023279312225
PARAF	903129486

Description

SCHIÖTZ TONOMETER - inclined scale

A masterpiece of precision mechanism, measures intraocular pressure according to Prof. Schiötz standards.
Vital parts are made of stainless steel.
High quality agate bearing for extremely long service life.
Set 3 weights included (5.5-7.5-10 g).
Scale specification 5.
In black deluxe case with velvet-look inserts.
Perfect reading of the scale with red pointer.

Made in Germany.