



Product Code	31245
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2733113
CND	V030208
EAN/UPC	8023279312454
UMDNS	14-255
PARAF	909314876
GMDN	32526

Description

ENT/NEUROLOGY TUNING FORKS

- Tuning fork a1-440 Hz - ENT - 90 g

All made of a special alloy and supplied with base.
 Multilingual manual and box: GB, FR, IT, ES, PT, GR.

Made in Germany.