



Product Code	31282
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2501265
CND	V030207
EAN/UPC	8023279312829
PARAF	963477310
GTIN	08023279312829
GMDN	11950

## Description

MONOFILAMENT TOOL - sensory evaluator

Plastic retractable monofilament with 3 spare monofilaments. Evaluates both diminishing and returning sensation (touch test). Examines levels of sensory function and changes in neuro status. Monofilament is calibrated at 10 g. Size: 13x2.5 cm + 4 cm monofilament

Diagnostic tool for detecting diabetic neuropathy.

Supplied in pouch with multilingual user manual (GB, FR, IT, ES, PT, DE, GR, Arabic) and 3 spare monofilaments