



| | |
|---------------|----------------|
| Product Code | 31289 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|---------|---------------|
| UK-REP | No |
| CH-REP | No |
| EAN/UPC | 8906200830602 |
| GMDN | 12293 |

Description

WORTH LED LIGHT TEST - with red/green glasses

The Worth 4 dot test, is a clinical test for suppression of either the right or left eye. Suppression occurs during binocular vision, when the brain does not process the information received from either of the eyes. This is a common adaptation to strabismus, amblyopia, and aniseikonia.

Battery not included, glasses included.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.