



Product Code	31290
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1171609
CND	Y210399
EAN/UPC	4580759310385
PARAF	977817648
GMDN	12293

## Description

**ISHIHARA COLOUR TEST**  
Book of 10 plates for illiterates.

Testing for colour deficiency requires high level of sensitivity and specificity.

This series of plates is provided to detect colour blindness in people who can not read, and colour deficiency, since young age. Plates are designed with simple geometrical shapes, since the examinees may not be able to understand complicated pieces.

Test gives quick and accurate assessment of colour vision deficiency of congenital origin.

Most cases are diagnosed with red-green deficiency which may be of two types:

- protan, which may be absolute (protanopia) or partial (protanomalìa)
- deutan, which may be absolute (deutanopia) or partial (deutanomalìa)