

ETHICON PERMA-HAND SILK SUTURES - gauge 4/0 needle 19 mm - braided - straight



Product Code	22313
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2470983
CND	H010202020301
GTIN	30705031088546
GMDN	13910

Description

PERMA-HAND - non absorbable sterile silk sutures, braided

- Ethicon code: 18507G
- Gauge: 4/0
- Metric: 1.5
- Colour: black
- Circle: Straight
- Needle length: 19 mm
- Needle: ST-4
- Suture length: 45 cm
- Minimum order: box of 12 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.