



Product Code	31314
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	7649996943064
---------	---------------

## Description

### LANG STEREOTEST® II-R

High and modern quality  
Improved optical quality of test cards thanks to modern digital and lenticular technology.

Homogeneity  
Increased homogeneity of the dot pattern and reduced single-eye recognition.

Rounded corners  
New surface lamination and laser cutting.

Correct application  
Schematic drawing and brief instructions on the back.

MDR 2017/745 compliance

LANG STEREOTEST® represents today's standard of random dot stereo tests without the need of 3D glasses, with very high sensitivity and specificity for stereopsis.  
The test combines random dots and a lenticular grid surface for quick screening of binocular vision in all groups of age, including very young children.  
While the test card is presented at reading distance, the patient is asked about what he sees and his eye movements are being observed.  
Recognition of figures is only possible if the person has full global stereopsis, while patients with amblyopia and (even small angle) strabismus fail the test.

LANG STEREOTEST® II-R has four stereo figures with disparities of 600, 400, 200 and 200 arcsec, of which one is also monocularly visible.

Multilingual instructions: GB, FR, IT, ES, PT, DE, PL, HU, NL, CZ, SE, DK, RO, HR, GR available on [www.gimaitaly.com](http://www.gimaitaly.com).

Made in Switzerland.