

ETHICON PERMA-HAND SILK SUTURES - gauge 4/0
needle 13 mm - braided



Product Code	22314
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2471082
CND	H010202020301
EAN/UPC	8023279223149
GMDN	13910

Description

PERMA-HAND - non absorbable sterile silk sutures, braided

- Ethicon code: 641G
- Gauge: 4/0
- Metric: 1.5
- Colour: black
- Circle: 3/8
- Needle lenght: 13 mm
- Needle: P-3 prime
- Suture lenght: 45 cm
- Minimum order: box of 12 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.