

XENON-HALOGEN VET DIAGNOSTIC SET



Product Code	31424
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279314243

Description

Includes: xenon halogen open otoscope head+handle, xenon-halogen opthalmoscope head, kit of 5 reusable plastic cannulas (31451), metal tongue depressor holder illuminator, nasal specula.

Rechargeable battery handle 2.5V

All sets have 2.5V handle for use with 2xC batteries or rechargeable batteries (34479)