



| | |
|---------------|----------------|
| Product Code | 31435 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2830656 |
| CND | Z12149006 |
| EAN/UPC | 8023279314359 |
| PARAF | 903146482 |
| GTIN | 08023279314359 |
| GMDN | 12849 |

Description

PARKER COLOUR OTOSCOPE - green
 Connection system: Bayonet
 Ear Funnels: Ø 2.5 - 3.5 - 4.5 mm

- Reliable and indestructible high quality instrument
- Invaluable for everyday ear nose and throat use
- Easily adjustable Neutral 2.5V vacuum light (made in Germany)
- Quality 3 times magnifying lens
- Available with Plug in and snap (bayonet locking) or knurled screw type otoscope heads
- Ideal for pneumatic tympanic membrane tests, can also be used for rhinoscopic examinations
- Supplied with 3 autoclavable ear funnels in a colour coordinated zip bag. All latex free
- Require 2 x C batteries (not included)

Multilingual box and manual: GB, FR, IT, ES, PT, GR, DE, Arabic.