



Product Code	31441
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2830656
CND	Z12149006
EAN/UPC	8023279314410
GTIN	08023279314410
GMDN	12849

Description

PARKER OTOSCOPE - black
Rechargeable handle
Connection system: Bayonet
Ear Funnel: Ø 2.5 - 3.5 - 4.5 mm

- Reliable and indestructible high quality instrument
- Invaluable for everyday ear nose and throat use
- Easily adjustable Neutral 2.5V vacuum light (made in Germany)
- Quality 3 times magnifying lens
- Available with Plug in and snap (bayonet locking) or knurled screw type otoscope heads
- Ideal for pneumatic tympanic membrane tests, can also be used for rhinoscopic examinations
- Supplied with 3 autoclavable ear funnels in a colour coordinated zip bag. All latex free
- Require 2 x C batteries (not included) and also works with rechargeable battery

Multilingual box and manual: GB, FR, IT, ES, PT, GR, DE, Arabic.