



Product Code	31442
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279314427

Description

PARKER VETERINARY OTOSCOPE - black
 Connection system: Screw
 Ear Funnels: Ø 4 - 5 - 7 mm

- Reliable and indestructible high quality instrument
- Invaluable for everyday ear nose and throat use
- Easily adjustable Neutral 2.5V vacuum light (made in Germany)
- Quality 3 times magnifying lens
- Available with Plug in and snap (bayonet locking) or knurled screw type otoscope heads
- Ideal for pneumatic tympanic membrane tests, can also be used for rhinoscopic examinations
- Supplied with 3 autoclavable ear funnels in a colour coordinated zip bag. All latex free
- Require 2 x C batteries (not included)

Multilingual box and manual: GB, FR, IT, ES, PT, GR, DE, Arabic.