

# KIT OF 5 REUSABLE PLASTIC SPECULA



Product Code	31451
Unit of sale	set of 5 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2654261
CND	Z12021007
EAN/UPC	8023279314519
GTIN	08023279314519
GMDN	33395

## Description

### AUTOCLAVABLE SPECULUM FOR OTOSCOPES - LATEX FREE

Kit of 5 autoclavable speculum Ø 2.5, 3.5, 4.5, 5.5, 9 mm

To be used with GIMA code: 31435-37, 31440-42, 31471, 31473, 31475, 31477, 31480  
Also compatible with: Heine (Non compatibile con Mini 3000), Kawe, Riester

- made from resilient thermoplastic material
- smooth matte finish surface ensure easy cleaning
- metallic bayonet lock for positive fit to otoscope
- autoclavable at 134°C