

KIT OF 3 REUSABLE METAL SPECULA + 2 METAL SPECULA



Product Code 31453

Unit of sale set of 5 pcs

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279314533

Description

AUTOCLAVABLE SPECULUM FOR OTOSCOPES - LATEX FREE

Kit of 3 metal veterinary cannulas (Ø 4, 5, 6 mm) + 2 slotted metal specula (Ø 7.5 and 9.4 mm)

To be used with GIMA code: 31422, 31424, 31425

2 slotted models with 60 mm opening and round edges, ideal for instrumentation, atraumatic insertion.
Matte inner surface to reduce risk of reflection. For use with Vet open otoscope heads