



Product Code	31467
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2681599
CND	Q030499
EAN/UPC	0014646063336
GMDN	65445

## Description

CERASPOON curettes: For viscous cerumen. For extensive cerumen buildup. For when a large scoop works better 4 mm wide tip.

Our four styles of curettes give you the most effective, economical and safe way to remove cerumen - whatever the size of ear or consistency of cerumen.

- Patient comfort: metal curettes are cold, hard and sharp. Bionix ear curettes with rounded plastic tips reduce trauma and risk of laceration. Latex free.
- Effectiveness: ear curettes can duplicate the effectiveness of metal curettes in almost any situation and their polypropylene composition makes them virtually unbreakable.
- Disposability: ear curettes don't require expensive and time-consuming sterilization and eliminate the risk of cross-contamination.