



Product Code	31469
Unit of sale	box of 48 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	214939
CND	L149004
EAN/UPC	0014646098772
GMDN	42449

## Description

DISPOSABLE NASAL SPECULUM - box of 48

Innovative plastic, disposable nasal speculum with handle: allows superior view of nasal cavity.

Helps extract foreign objects from nasal cavity.  
Facilitates packaging of cauterizing of bleeding noses.  
Latex free.

Made in USA.

The use of stainless steel speculum can be a chilling experience for the patient. Now you can combine patient comfort with effective treatment. Smooth, warm tips make disposable Nasal Speculum more comfortable and less alarming to the patient.

Total blade depth: 33 mm  
Total length: 13.2 cm