



Product Code	31471
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2830688
CND	Z121402
EAN/UPC	8023279314717
PARAF	903626036
GTIN	08023279314717
GMDN	12849

Description

PARKER DIAGNOSTIC SET

Lightweight, rigid case with one Parker otoscope 2.5 V with rechargeable handle, 3 autoclavable ear speculum (Ø 2.5, 3.5 and 4.5 mm), 2 laryngeal mirrors (size 3 and 4), bent arm throat lamp, tongue spatula, adjustable nasal speculum.