



Product Code	31472
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2681599
CND	Q030499

Description

VERSALOOP® - purple - box of 50 pcs

VERSALOOP®: for gentle precision Smooth and flexible tip that gives the user exceptional control while offering the patient more comfort. 3 mm tip

Bionix Safe Ear Curettes™ allow healthcare providers to safely and gently remove cerumen with styles to match different types of wax and patient.

Patient comfort: Bionix Safe Ear Curettes are flexible, to reduce injury and pain to the ear canal, making them safer than hard and rigid stainless-steel currettes.

Effectiveness: Safe Ear Curettes can be more effective than stainless steel currettes. They are made of smooth, flexible polypropylene designed to bend when there's an obstacle.

Disposability: single-use currettes reduce the risk of crosscontamination and eliminate sterilization costs and hassles.

Made in the USA.