



Product Code	31475
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2830656
CND	Z12149006
EAN/UPC	8023279314755
PARAF	903146684
GMDN	12849

Description

Black model with halogen 2.5 V bulb and bayonet connection system.

- Reliable and indestructible high quality instrument
- Invaluable for everyday ear nose and throat use
- Easily adjustable Neutral 2.5V vacuum light (made in Germany)
- Quality 3 times magnifying lens
- Ideal for pneumatic tympanic membrane tests, can also be used for rhinoscopic examinations

Supplied with 3 autoclavable ear funnels in a soft bag. Requires 2xC batteries (not included) or rechargeable battery.

Multilingual box and manual: GB, FR, IT, ES, PT, GR, DE, Arabic.