

INFLATION BULB FOR OTOSCOPES - grey



| | |
|---------------|----------------|
| Product Code | 31483 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2324878 |
| CND | Q0399 |
| EAN/UPC | 8023279314830 |

Description

Pneumatic PVC bulb for Oscopes

To be used with GIMA code: 31435/36/37/40/42/71/73/75/77-31520/21/22, 31535
Also compatible with: Heine, Kawe, Riester

Weight 60g
Total length 54 cm, tube 46 cm