

MINI EAR SPECULUM diam. 2.5 mm - grey - in dispenser



Product Code	31485
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2611667
CND	Q030499
EAN/UPC	8023279314854
UMDNS	10571
PARAF	951157609
GTIN	08023279314854
GMDN	33395

Description

GIMA DISPOSABLE SPECULUM FOR MINI OTOSCOPES - LATEX FREE

- Disposable Mini ear-speculum Ø 2.5 mm - grey

Compatibility

- Direct fit:

Gima Gimalux, Sigma, Sigma C, Visio 5, Visio 2000
Heine mini 2000, 3000, K180/Beta 200/Beta 400
Riester Ri-Mini, Ri-scope L1-L2, E-scope, Pen-scope
Kawe (Picolight, Combilight F030, Eurolight F030)

Holtex
Colson
Luxamed

- Fit with connector 31493:

Gima: 31435-37, 31440-42, 31471, 31473, 31475, 31477, 31480, 31493 (Parker line)