



Product Code	31494
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	151893
CND	L160301
EAN/UPC	8023279314946
PARAF	903146835
GTIN	08023279314946
GMDN	33395

Description

Alligator pediatric connector green colour, supplied with 100 green speculum Ø 2.5 mm

To be used with GIMA codes:

31435, 31436, 31437, 31440, 31441, 31442, 31471, 31473, 31475, 31477, 31480