



Product Code	31495
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2611667
CND	Q030499
EAN/UPC	8023279314953
UMDNS	10571
PARAF	903146759
GTIN	08023279314953
GMDN	33395

## Description

### GIMA DISPOSABLE SPECULUM FOR MINI OTOSCOPES - LATEX FREE

- Disposable Mini ear-speculum mixed colours and Ø
- Includes 500 ear speculum Ø 4 mm and 500 Ø 2.5 mm (200 black, 100 red, blue, green each size)

#### Compatibility

- Direct fit:

Gima Gimalux, Sigma, Sigma C, Visio 5, Visio 2000  
Heine mini 2000, 3000, K180/Beta 200/Beta 400  
Riester Ri-Mini, Ri-scope L1-L2, E-scope, Pen-scope  
Kawe (Picolight, Combilight F030, Eurolight F030)

Holtex  
Colson  
Luxamed

- Fit with connector 31493:

Gima: 31435-37, 31440-42, 31471, 31473, 31475, 31477, 31480, 31493 (Parker line)