

# HEINE MINI EAR SPECULUM - diameter 2.5 mm - grey



Product Code	31502
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	85569
CND	Z12149085
EAN/UPC	4053755114160
UMDNS	13-662

## Description

HEINE DISPOSABLE MINI EAR SPECULUM  
Compatible with Heine Mini 2000/3000, K180, Beta 200, Beta 400

Disposable Allspec® spec. Ø 2.5 mm - grey - dispenser of 250  
HEINE code: B-000.11.151