

**EAR SPECULUM diam. 4 mm for WelchAllyn - disposable**  
**- black**



|               |                |
|---------------|----------------|
| Product Code  | 31509          |
| Unit of sale  | box of 250 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2078999        |
| CND        | Q030499        |
| EAN/UPC    | 3701126705318  |
| GTIN       | 03701126705318 |

**Description**

Disposable ear-speculum Ø 4 mm - black

To be used with:  
Welch Allyn Pocket-Pro and Pocket scope Series 200, 250

Multilingual box: GB, FR, IT, ES, PT, DE, GR, DK, Arabic.