

# SIGMA F.O. OTOSCOPE 2.5V - Xenon-halogen light - pouch - violet



Product Code	31511
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2831178
CND	Z12149006
EAN/UPC	8023279315110
GTIN	08023279315110
GMDN	12849

## Description

- SIGMA F.O. otoscope 2.5V - standard handle
- Violet
- Xenon-Halogen
- Packing: Pouch

Supplied with 14 disposable specula (7 with Ø 2.5 mm, 7 with Ø 4 mm).

High quality pocket diagnostic instruments designed for professional use in ENT and ophthalmoscopy. Light weight and manufactured from impact resistant re-enforced plastic.

Manufactured with high quality materials and top-notch workmanship ensures excellent reliability in every examination.

- distal fibre optic illumination with super bright Xenon/Halogen light
- swivel type window with 3X magnification lens for sharp images and minimal reflection
- insufflation port for pneumatic testing of tympanic mobility

Connector port (code 31532) and insufflation bulb (code 31500) should be ordered separately.