

HEINE ALLSPEC MINI EAR SPECULUM - diameter 4 mm - grey



Product Code	31516
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	85569
CND	Z12149085
EAN/UPC	4053755113897
GMDN	34897

Description

HEINE DISPOSABLE MINI EAR SPECULUM
Compatible with Heine Mini 2000/3000, K180, Beta 200, Beta 400

Disposable Allspec® spec. Ø 4 mm - grey - box of 1000
HEINE code: B-000.11.127