



Product Code	31522
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2835632
CND	Z12149006
EAN/UPC	8023279315226
PARAF	903146898
GTIN	08023279315226
GMDN	12849

## Description

GIMALUX F.O. OTOSCOPE - green

High professional plastic otoscopes with extra bright halogen light provided by fibre optic system. Quality 3 times magnifying movable lens. Output for pneumatic membrane test. Supplied with 14 single use mini speculums (7 Ø 2.5 and 7 Ø 4 mm) in a black nylon bag. All latex free.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, GR, RO, HU, Arabic.