

SIGMA F.O. LED OTO-OPHTHALMOSCOPE SET with 2 handles - case



Product Code	31534
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2734775
CND	L1699
EAN/UPC	8023279315349
GMDN	45084

Description

SIGMA F.O. LED OTO-OPHTHALMOSCOPE DELUXE SET 2.5 V - case with 2 handles

Oto-ophthlmo sets are supplied with oto and ophtalmo head, 2 handles and 14 disposable speculum (7x2.5 and 7x4 mm).