



Product Code	31535
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2831188
CND	Z12149006
EAN/UPC	8023279315356
UMDNS	12849
GMDN	12849

## Description

VISIO 2000 F.O. XENON OTOSCOPE - 3.5 V

F.O. otoscope head, battery handle with bayonet lock, 1 set of 5 reusable eartips (2.5, 3.0, 4.0, 5.0 & 10.0 mm), tube of 14 disposable eartips (2.5 & 4.0 mm), spare bulb, in a rigid case.

Professional diagnostic sets with fiber optic illumination, Xenon light 3.5 V, rheostat to adjust the intensity of light.  
Made of resistant re-inforced plastic and durable easy locking system.  
Supplied in deluxe case, should be used with charging station 31542 and Li-ion rechargeable battery 31541.

## Technical Specifications

Xenon otoscope

- Swivel type window with 3X magnification for sharp images and minimal reflection
- Distal F.O. illumination with super bright LED light
- Durable easy locking system
- Rheostat