



Product Code	31536
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2831204
CND	Z12120114
EAN/UPC	8023279315363
UMDNS	12849
GTIN	08023279315363
GMDN	45084

Description

VISIO 2000 F.O. XENON OPHTHALMOSCOPE - 3.5 V

F.O. ophthalmo head, battery handle with bayonet lock, spare bulb, in a rigid case.

Professional diagnostic sets with fiber optic illumination, Xenon light 3.5 V, rheostat to adjust the intensity of light. Made of resistant re-inforced plastic and durable easy locking system. Supplied in deluxe case, should be used with charging station 31542 and Li-ion rechargeable battery 31541.

Technical Specifications

Xenon ophthalmoscope

- 6 apertures to select from:
 - Micro spot aperture
 - Small aperture
 - Large aperture
 - Fixation aperture
 - Slit
 - Red-free filter
- Enhanced Xenon super bright light
- Durable easy locking system
- Rheostat