

VISIO 2000 F.O.XENON OTO-OPHTHALMOSCOPE PLUG IN - 3.5 V



Product Code	31537
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2734804
CND	L1699
EAN/UPC	8023279315370
GMDN	45084

Description

VISIO 2000 F.O. XENON OTO-OPHTHALMOSCOPE PLUG-IN SET - 3.5 V - 2 pin plug

F.O. otoscope head, F.O. ophthalmo head, battery handle with bayonet lock and plug-in, 1 set of 5 reusable eartips (2.5, 3.0, 4.0, 5.0 & 10.0 mm), tube of 14 disposable eartips (2.5 & 4.0 mm), spare oto and ophthalmo bulbs, rechargeable Li-ion battery 31547, in a rigid case.

Professional diagnostic sets with fiber optic illumination, Xenon light 3.5 V, rheostat to adjust the intensity of light.

Made of resistant re-inforced plastic and durable easy locking system.

Supplied in deluxe case, should be used with charging station 31542 and Li-ion rechargeable battery 31541.