



|               |                |
|---------------|----------------|
| Product Code  | 31539          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2734804       |
| CND        | L1699         |
| EAN/UPC    | 8023279315394 |
| UMDNS      | 12849 + 46788 |
| GMDN       | 45084         |

## Description

### VISIO 2000 F.O. XENON DIAGNOSTIC SET - 3.5 V

F.O. otoscope head, F.O. ophthalmo head, battery handle with bayonet lock, nasal illuminator with F.O. light, curved illuminator with bulb, plain mirror 20.0 mm, plastic tongue depressor, 1 set of 5 reusable eartips (2.5, 3.0, 4.0, 5.0 & 10.0 mm) tube of 14 disposable eartips (2.5 & 4.0 mm), spare oto and ophthalmo bulbs, in a rigid case.

Professional diagnostic sets with fiber optic illumination, Xenon light 3.5 V, rheostat to adjust the intensity of light. Made of resistant re-inforced plastic and durable easy locking system. Supplied in deluxe case, should be used with charging station 31542 and Li-ion rechargeable battery 31541.

## Technical Specifications

### Xenon otoscope

- Swivel type window with 3X magnification for sharp images and minimal reflection
- Distal F.O. illumination with super bright LED light
- Durable easy locking system
- Rheostat

### Xenon ophthalmoscope

- 6 apertures to select from:
  - Micro spot aperture
  - Small aperture
  - Large aperture
  - Fixation aperture
  - Slit
  - Red-free filter
- Enhanced Xenon super bright light
- Durable easy locking system
- Rheostat