

E-SCOPE OTOSCOPE - halogen 2.7V - black in pouch



Product Code	31548
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	1830774
CND	Z12149006
EAN/UPC	4045396176898

Description

E-SCOPE® OTOSCOPE - HALOGEN 2.7V - BLACK IN POUCH

Halogen illumination otoscope.
Sealing precision optics pivoting to both sides with triple magnification for superior view.
Suitable for performance of pneumatic tests (supplied without connector and insufflator).