

# E-SCOPE F.O. OTOSCOPE - LED 3.7V - white in rigid case



Product Code	31552
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396176928
---------	---------------

## Description

### E-SCOPE® F.O. OTOSCOPE - LED 3.7V - WHITE IN RIGID CASE

F.O. LED illumination otoscope.

Sealing precision optics pivoting to both sides with triple magnification for superior view.

Suitable for performance of pneumatic tests (supplied without connector and insufflator).

#### INNOVATIVE LED 3.7 V ILLUMINATION :

- real colour 5,500K white light

- long life LED bulbs: up to 20,000 hours

- low battery consumption: up to 50 hours operation through IPC (Integrated Power Converter). It converts the voltage of AA batteries into the optimum voltage required for LED