

E-SCOPE OPHTHALMOSCOPE - Xenon 2.5V - black in rigid case



Product Code	31556
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396176997
---------	---------------

Description

E-SCOPE® OPHTHALMOSCOPE - XENON 2.5 V - BLACK - IN RIGID CASE

- with 2.5V Xenon
- dioptre disc with 18 corrective lenses ($\pm 1, 2, 3, 4, 6, 8, 10, 15, 20$)
- easy-to-use aperture wheel with six apertures (fixation star, large circle, small circle, red free filter, blue filter, semi circle)
- high performance optics with aspherical condenser
- spectacle protection
- parallel optical path