

# E-SCOPE DIAGNOSTIC KIT - LED 3.7V - black in rigid case



Product Code	31564
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

## Description

### E-SCOPE® DIAGNOSTIC KIT - BLACK

Oto-ophthalmoscope sets, LED 3.7V, in rigid case.

#### • Otoscope

Sealing precision optics pivoting to both sides with triple magnification for superior view. Suitable for performance of pneumatic tests (supplied without connector and insufflator).

#### • Ophthalmoscope

- innovative 3.7V LED illumination for higher-contrast and more efficient diagnosis
- dioptre disc with 18 corrective lenses ( $\pm 1, 2, 3, 4, 6, 8, 10, 15, 20$ )
- easy-to-use aperture wheel with six apertures (fixation star, large circle, small circle, red free filter, blue filter, semi circle)
- high performance optics with aspherical condenser
- spectacle protection
- parallel optical path

#### INNOVATIVE LED 3.7 V ILLUMINATION :

- real colour 5,500K white light
- long life LED bulbs: up to 20,000 hours
- low battery consumption: up to 50 hours operation through IPC (Integrated Power Converter). It converts the voltage of AA batteries into the optimum voltage required for LED