

LUXAMED AURIS LED 2.5V OTOSCOPE - black



Product Code	31570
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2219532
CND	Z12149006
EAN/UPC	8023279315707

Description

LUXAMED AURIS LED 2.5V OTOSCOPE - BLACK

Innovative otoscopes that combine best handicraft production with selected materials with the most modern technology. Two-components handle (aluminium/fibre-glass reinforced) with electronic ON/OFF switch and dimmer.

- Features:
- LED-Ring illumination with 6 LEDs (concentric), no fibre optic. The functionality is not affected (patent DE202013104278U1) if a single LED fails.
 - automatic switchoff after 3 minutes
 - approx. 100,000 hours life expectancy of LEDs
 - approx. 10,000 Lux intensity of illumination
 - approx. 4,000 K colour temperature
 - pivoting 3-fold magnifying glass with locking function
 - insufflation port for pneumatic test.

Supplied with 2xAAA alkaline batteries, 20 disposable ear funnels (10x2.5 mm and 10x4.0 mm) and manual (GB, IT, DE).

Made in Germany.