



Product Code	31572
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2219532
CND	Z12149006
EAN/UPC	8023279315721

Description

LUXAMED AURIS LED 3.7V USB OTOSCOPE - BLACK

USB port: recharger with USB plug - 3 adapters (EU, UK, US)

Innovative otoscopes that combine best handicraft production with selected materials with the most modern technology. Two-components handle (aluminium/fibre-glass reinforced) with electronic ON/OFF switch and dimmer.

Features:

- LED-Ring illumination with 6 LEDs (concentric), no fibre optic. The functionality is not affected (patent DE202013104278U1) if a single LED fails.
- automatic switchoff after 3 minutes
- approx. 100,000 hours life expectancy of LEDs
- approx. 10,000 Lux intensity of illumination
- approx. 4,000 K colour temperature
- pivoting 3-fold magnifying glass with locking function
- insufflation port for pneumatic test.

Supplied with lithium-ion rechargeable battery 3.7V, 20 disposable ear funnels (10x2.5 mm and 10x4.0 mm) and manual (GB, IT, DE).

Made in Germany.